

SwabPro™

Every Drop Matters



How SwabPro™ Works

1

ATTACH

Attach to an
ENFit® syringe

2

ADMINISTER

Controlled, direct
oral delivery.



IDEAL FOR

Rx

Buccal/
Oral Delivery



Milk Drops



Colostrum
Oral Care



Glucose
Delivery

CONTROLLED • GENTLE • VERSATILE

See how swabpro compares to conventional methods.

SwabPro™ at a Glance

CATEGORY	SwabPro™ 	Gloved Finger 	Syringe Only 	Cotton Swab Dip 
 Dosing Accuracy	✓	✗	✓	✗
 Application Control	✓	✗	✗	✗
 Consistency Between Caregivers	✓	✗	✗	✗
 Hygiene & Safety	✓	✗	✓	✗
 Gentleness to Oral Tissue	✓	✓	✗	✓
 Infant-Led Interaction	✓	✓	✗	✓
 Non-Shedding Foam Tip	✓		✗	✗



= Meets standard/practice



= Does not meet


Every Drop Matters



One Device



Multiple Oral Therapies



Consistent Bedside Practice



Gentle Controlled Delivery

Ordering Information

Product Code	Description	Packaging	Case Quantity
SP&NS	SwabPro™ Oral-Therapy Applicator	30 each per box	600 per case (20boxes of 30)

SwabPro™
Every Drop Matters

